
✅ Chameleon KVM: End-to-End MLOps Pipeline Workflow

This guide walks through setting up a scalable MLOps system for chest X-ray diagnosis using Docker containers, persistent storage, and a feedback loop, deployed on a manually configured Chameleon KVM with assigned public IP.

1. 🗝️ Access KVM Instance

Start by SSH-ing into your KVM instance:

```
ssh -i ~/.ssh/id_rsa_chameleon_project44 cc@<KVM_IP>
```

2. 📁 Clone the Project Repository

```
git clone https://github.com/Shrey12202/MedAI-Scalable-Diagnosis-with-Machine-Learning.git
```

3. 🐳 Launch Core Infrastructure Containers

Start services: **MLflow**, **Grafana**, **Prometheus**, **MinIO**, and **Postgres**.

```
cd MedAI-Scalable-Diagnosis-with-Machine-Learning/CheXspert\MLOps/kvm\ container
```

```
docker compose -f docker-compose-mlflow-minio-grafana-prometheus-postgres up -d
```

Ports Used:

Service	Port
MLflow	8000
Grafana	3000
Prometheus	8001
MinIO	9000
Postgres	9001

4. 🚀 Launch Inference & Feedback Services

Start **FastAPI**, **Label Studio**, and the **feedback transfer** component:

```
cd MedAI-Scalable-Diagnosis-with-Machine-Learning/CheXspert\MLOps/serving/docker
```

```
docker compose -f eval-loop-chi/docker/docker-compose-production.yaml up -d
```

```
docker compose -f eval-loop-chi-main/docker/docker-compose-labelstudio.yaml up -d
```

```
docker build -t user-feedback-transfer .
```

```
docker run --rm user-feedback-transfer
```

Service	Port
FastAPI	5000
Label Studio	8080

5. 🗄️ Mount Persistent Storage (CHI@TACC)

On a node with access to node-persist, run:

```
curl https://rclone.org/install.sh | sudo bash
```

```
sudo sed -i '/^#user_allow_other/s/^#//' /etc/fuse.conf
```

```
mkdir -p ~/.config/rclone
```

```
nano ~/.config/rclone/rclone.conf
```

Paste into rclone.conf:

```
[chi_tacc]
```

```
type = swift
```

```
user_id = sdv2034
```

```
application_credential_id = 292824fe1f334db69bc98130063de7f1
```

```
application_credential_secret = <your-secret>
```

```
auth = https://chi.tacc.chameleoncloud.org:5000/v3
```

```
region = CHI@TACC
```

Then mount the object store:

```
sudo mkdir -p /mnt/object
```

```
sudo chown -R cc /mnt/object
```

```
sudo chgrp -R cc /mnt/object
```

```
rclone mount chi_tacc:object-persist-project44-1 /mnt/object \  
--read-only --allow-other --vfs-cache-mode=full \  
--dir-cache-time=72h --swift-fetch-until-empty-page --daemon
```

6. 🧠 Model Training

Prepare Docker and run the training pipeline:

```
sudo apt update && sudo apt install docker.io
```

```
sudo systemctl start docker && sudo systemctl enable docker
```

```
cd MedAI-Scalable-Diagnosis-with-Machine-Learning/CheXspert\ MLOps/train
```

Optional: Add user to Docker group

```
sudo usermod -aG docker $USER && newgrp docker
```

```
docker build -t chexpert-train .
```

```
docker run --shm-size=16g \  
-v $(pwd):/workspace \  
-v /mnt/object:/mnt/object \  
chexpert-train
```

```
-v /mnt/persistent:/mnt/persistent \
chexpert-train
```

7. Monitoring, Inference & Feedback

- **MLflow UI (metrics, models):** <http://129.114.27.181:8000>
 - **MinIO (artifacts):** <http://129.114.27.181:9000>
 - **FastAPI (inference):** <http://129.114.27.181:5000>
 - **Label Studio (review/correct):** <http://129.114.27.181:8080/projects/>
 - **Grafana (data drift & metrics):** <http://129.114.27.181:3000>
-

8. Feedback Transfer & Retraining Loop

After running inference:

```
cd MedAI-Scalable-Diagnosis-with-Machine-Learning/CheXspert\ MLOps/serving/workspace
```

```
docker build -t user-feedback-transfer .
```

```
docker run --rm user-feedback-transfer
```

This moves flagged predictions into Label Studio for manual correction. Once corrected, retraining can begin again from step 6, using:

- Original data from /mnt/object
 - Simulated feedback data stored in MinIO
-